




# Ultraviolet-Visible Spectroscopy









**Figure 13-1** Reflection and scattering losses.

# NATURE OF ELECTRONIC EXCITATIONS

Absorption: A transition from a **lower level** to a **higher level** with **transfer of energy** from the **radiation field** to an **absorber, atom, molecule, or solid**.


Emission: A transition from a **higher level** to a **lower level** with **transfer of energy** from the **emitter** to the **radiation field**.






## ❑ The Spectroscopic Process

1. In **UV spectroscopy**, the sample is **irradiated** with the broad spectrum of the **UV radiation**
2. If a particular **electronic transition** matches the **energy** of a certain **band of UV**, it will be **absorbed**
3. The remaining **UV light** passes through the sample and is **observed**
4. From this residual radiation a **spectrum** is obtained with “**gaps**” at these **discrete energies** – this is called an **absorption spectrum**.



## □ Observed electronic transitions

From the **molecular orbital diagram**, there are several possible **electronic transitions** that can occur, each of a **different relative energy**



$\sigma \rightarrow \sigma^*$  **alkanes**

$\sigma \rightarrow \pi^*$  **carbonyls**

$\pi \rightarrow \pi^*$  **unsaturated compounds**

$n \rightarrow \sigma^*$  **O, N, S, halogens**

$n \rightarrow \pi^*$  **carbonyls**

# Transitions

## $\sigma \rightarrow \sigma^*$

UV photon required, high energy

Methane at 125 nm ( $\text{CH}_4$ )

Ethane at 135 nm ( $\text{C}_2\text{H}_6$ )

## $n \rightarrow \sigma^*$

Saturated compounds with unshared  $e^-$

Absorption between 150 nm to 250 nm

$\epsilon$  between 100 and 3000  $\text{L cm}^{-1} \text{ mol}^{-1}$

Shifts to shorter wavelengths with polar solvents

Minimum accessibility

Halogens, N, O, S

## $n \rightarrow \pi^*$ , $\pi \rightarrow \pi^*$

Organic compounds, wavelengths 200 to 700 nm

Requires unsaturated groups

$n \rightarrow \pi^*$  low  $\epsilon$  (10 to 100)

Shorter wavelengths

$\pi \rightarrow \pi^*$  higher  $\epsilon$  (1000 to 10000)

# INSTRUMENTATION

**Spectrometer:** An instrument used for measuring transmittance or absorbance of a sample as function of the wavelength of the electromagnetic radiation



## Components of a spectrophotometer

**Source:** It generate electromagnetic radiation

**Dispersion device:** It select the wavelength required from the broad band of radiation source

**Sample area:** Where the sample is kept

**Detector:** One or more detectors that measure the intensity of the radiation



# Radiation Sources


Two sources are required to scan the entire UV-Visible band:

- **Deuterium arc** lamp– covers the **UV – 200-330**
- **Tungsten-halogen** lamp – covers **330-700** (yields good intensity over part of the UV spectrum and over the entire visible range)
- An alternate light source: **Xenon** lamp
  - Pros:** Yields a good continuum over the entire **UV** and **visible** regions.
  - Cons:** **High noise** from currently available **Xenon lamps** compared to **deuterium** or **tungsten** lamps

□ **Monochromator:** consists of an **entrance** slit, a **dispersion** device, and an **exit** slit.



# Dispersion devices

- ❖ Dispersion devices cause different wavelengths of light to be dispersed at different angles.
- ❖ When combined with an appropriate exit slit, these devices can be used to select a particular wavelength of light from a continuous source.
- ❖ Two types of commonly used dispersion devices:
  1. Prisms
  2. Holographic gratings



# Sample array

1. Sample cells can be made of **plastic, glass or quartz**
2. **Glass** absorbs strongly below **320 nm**
3. The cells **lowest in cost** are made of **plastic**, usually an **acrylic**. These cells **are not resistant** to all solvents and **absorb strongly** below **300 nm**
4. Only **quartz** is transparent in the full **200-700 nm** range; **plastic** and **glass** are only suitable for **visible regions spectra**




# DETECTORS


❖ A **detector** converts a **light signal** into an **electrical signal**.

❖ **Spectrophotometers** normally contain

**Photomultiplier tube detector**  
or a  
**Photodiode detector.**




**Photomultiplier tube**



**Photodiode**

# SPECTROMETER DESIGN




## Conventional Spectrometer

Polychromatic source of light falls on an **entrance slit** which transmits a **narrow band** of light.

The light then pass through the **sample** to a detector.

The detector measures the **absorbance** of the sample **by comparing the light that reaches the detector from the sample and the blank (only solvent)**



## Diode array Spectrometer

Polychromatic source of light falls on a the sample, the transmitted radiation pass through an **entrance slit** of the dispersion device.

The detector measures the **absorbance** of the sample **by comparing the light that reaches the detector from the sample and the blank (only solvent)**